Bilirrubina Total

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ENGLISH

INTENDED USE

Kit intended to determine total bilirubin in serum and plasma (heparin and EDTA) samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: DCA (Diazotized Dichloroaniline).

In acidic pH, direct and indirect bilirubin are coupled with diazotized dichloroaniline, being solubilized and producing a colored complex, called azobilirubin, that can be measured at 550 nm (505-580 nm). The color intensity is proportional to the concentration of total bilirubin in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (heparin and EDTA).

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma*	4 to 8 °C	7 days
	-20 °C	6 months

*Necessary to keep out of the light when stored for longer than 8 hours.

PRODUCT DESCRIPTION

R 1 Surfactants, stabilizers, catalyst and detergent

R 2 Diazotized dichloroaniline < 2,0 mmol/L and stabilizer.
pH < 5

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal H	13.002,00
Normal Control - Qualinorm	13.003,00
Pathological Control - Quantialt	13.004,00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 550 nm (505 - 580 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

Blank	Sample Blank	Sample	Calibrator Blank	Calibrator
---	---	---	50 µL	50 µL
---	50 µL	50 µL	---	---
50 µL	---	---	---	---
R1	1.0 mL	1.0 mL	1.0 mL	1.0 mL
	Homogenize and incubate the test tubes for 1 minute at 37 °C			
R2	250 µL	---	250 µL	---

2. Homogenize and incubate at 37 °C for 4 minutes.

3. Measure the Blank's, Sample Blank's, Sample's and Standard's absorbance against the purified water at 550 nm (505-580 nm). The final reaction is stable for 30 minutes.

B) CALCULATIONS

Total Bilirubin (mg/dL) = $\frac{\text{Sample's } \Delta A - \text{Blank}}{\text{Calibrator's } \Delta A - \text{Blank}}$

In which:

ΔA = Sample's or Calibrator's Absorbance – Sample or Calibrator Blank's Absorbance

0.91	0.02	1.80	0.91	0.02	1.80
6.02	0.10	1.70	5.98	0.12	2.00
10.80	0.10	0.90	10.80	0.08	0.80

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RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Newborns: 1 day	Up to 5.1 mg/dL
Newborns: 1 to 2 days	Up to 7.2 mg/dL
Newborns: 3 to 5 days	Up to 10.3 mg/dL
Children, Adolescents and Adults	Up to 1.2 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicalda.com.br.

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