Tira de Urina

Urine Test Strip / Tira de Orina
Ref. 30.009.00

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ENGLISH

INTENDED USE

Kit for qualitative and semi-quantitative determination of glucose, bilirubin, ketone, density, blood, pH, protein, urobilinogen, nitrite and leukocytes in the urine. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 15 to 30 °C.
- The product is stable until the expiration date printed on the label, as long as the recommended storage conditions are followed.
- Once the bottle is opened, the strips are stable for 3 months, provided they are kept in the original well capped bottle and stored at a temperature of 15 to 30 ° C.
- Ready-to-use strip.

WORKING PRINCIPLE

Leukocytes (LEU): the test reveals the presence of granulocyte esterase. The esterase is reacted with the amino pyrazole acid ester releasing hydroxypyrazole. Pyrazole reacts with diazonium salt to produce a violet color. The intensity of color developed is proportional to the number of leukocytes present in the urine sample.

Nitrite (NIT): specific test for nitrites, which are formed by the reduction of nitrates due to the action of the reductases produced by Gram negative bacteria. In an acidified medium, the nitrite in the urine reacts with p-arsanilic acid to form a diazonic compound. This is coupled with 1N- (1-naphthyl) ethylenediamine to give a pink color.

Urobilinogen (URO): Ehrlich's modified colorimetric reaction between p-diethylaminobenzaldehyde and urobilinogenenic acid in strongly acid medium that produces a pink color.

Protein (PRO): in the presence of protein the blue indicator of tetrabromophenol changes color, varying from yellow to green.

pH: double methyl red / bromothymol blue indicator system. The results vary from orange to bluish-green.

Blood (BLO): the activity of hemoglobin pseudoperoxidase catalyzes the reaction of dihydroperoxydiiodisopropylbenzene and 3,3', 5,5'-tetramethylbenzidine. A uniform color ranging from yellow to blue indicates the presence of myoglobin, hemoglobin or hemolyzed erythrocytes. Scattered blue-green spots or compact spots indicate intact erythrocytes as scale of results.

Density (SG): Urine density is related to ionic concentration, based on the apparent change of pKa. The resulting color ranges from dark blue green to yellowish green.

Ketones (KET): Ketones react with nitroprusside and acetoacetic acid that produce a change in color, ranging from light pink to purple.

Bilirubin (BIL): Bilirubin binds with dichloroanilina diazotized in a strongly acidic medium. The variation in bilirubin levels will produce a color with intensity proportional to its concentration in the urine.

Glucose (GLU): Glucose is oxidized by glucose oxidase forming gluconic acid and hydrogen peroxide. Hydrogen peroxide reacts with the chromogen potassium iodide in the presence of peroxidase. The amount of oxidized chromogen determines the color produced, ranging from green to brown.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: urine

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation: The test should be performed on a fresh urine sample, with no preservative added and kept at room temperature. When analysis cannot be performed within two hours of collection, the sample should be refrigerated and protected from light, however refrigeration may increase specific gravity and precipitate amorphous phosphates and urates. Refrigeration conditions should be defined by the laboratory, taking into account local characteristics. The sample should never be frozen. Wait for urine to reach room temperature before analysis.

PRODUCT DESCRIPTION

Reagent Area	Composition
Leukocytes (LEU)	Amino pyrazole acid ester derivative; diazonium salt; buffer solution; non-reactive ingredients.
Nitrite (NIT)	p-arsanilic acid; N- (1-naphthyl) ethylenediamine; non-reactive ingredients.
Urobilinogen (URO)	p-diethylaminobenzaldehyde; buffer solution and non-reactive ingredients.
Protein (PRO)	Tetrabromophenol blue; buffer solution and non-reactive ingredients.
pH	Sodium methyl red salt; bromothymol blue; non-reactive ingredients.
Blood (BLO)	3,3', 5,5'-tetramethylbenzidine (TMB); diisopropylbenzene dihydroperoxide. Buffer solution and non-reactive ingredients.
Density (SG)	Blue bromothymol indicator; buffer solution and non-reactive ingredients; poly (methyl vinyl ether / anhydride maleic) sodium hydroxide.
Ketones (KET)	Sodium nitroprusside; solution.
Bilirubin (BIL)	2,4-dichloroaniline diazonium salt; buffer solution and non-reactive ingredients.
Glucose (GLU)	Glucose oxidase; peroxidase; potassium iodide; buffer solution; non-reactive ingredients.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use control in the reference range or at the decision and control level with value in another range of clinical significance

NECESSARY EQUIPMENT FOR TESTING


- Container for sample collection;
- Clock or stopwatch.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Homogenize the urine.
2. Remove the strips and close the tube immediately.
3. Immerse test strip in urine for approximately 2 seconds.
4. Remove the strip by sliding it over the edge of the sample tube to remove excess urine.
5. Hold the strip in the horizontal position and put it in contact with absorbent paper to avoid mixing chemical reagents from adjacent reaction areas.
6. Wait for the time specified below and read the results:

Glucose and bilirubin	30 seconds
Ketones	40 seconds
Density	45 seconds
Blood, pH, protein, urobilinogen and nitrite	60 seconds
Leukocytes	120 seconds

B) ANALYSIS OF RESULTS

The results are obtained by visual comparison with the color scales printed on the vial label.

C) INTERPRETATION

Leukocytes (LEU): Normal urine samples present negative results.

Nitrite (NIT): Nitrite is not present in normal urine.

Urobilinogen (URO): Normal urine contains 3.5-17 µmol / L of urobilinogen.

Protein (PRO): Normal urine shows a negative reaction.

pH: 5.0-7.0.

Blood (BLO): Normal urine samples present negative results.

Density (SG): Between 1,000 and 1,030.

Ketones (KET): Normal urine samples present negative results.

Bilirubin (BIL): Normal urine samples present negative results.

Glucose (GLU): Normal urine shows negative reaction.

PERFORMANCE CHARACTERISTICS

Operating Range	
Leukocytes (LEU)	Negative (-) at 500 Leu/µL (+++).
Nitrite (NIT)	Negative and positive.
Urobilinogen (URO)	0.2 to 12 mg/dL (3.5 to 200 µmol/L).
Protein (PRO)	Negative at 2000 mg/dL (20 g/L).
pH	5.0 to 9.0.
Blood (BLO)	Negative at (+++); 5-10 to 50 Ery/µL.
Density (SG)	1.000 to 1.030.
Ketones (KET)	Negative at 160 mg/dL (16 mmol/L).
Bilirubin (BIL)	Negative at 4.0 mg/dL (70 µmol/L).
Glucose (GLU)	Negative at ≥ 2000 mg/dL (110 mmol/L).

Sensitivity			
Leukocytes	Nitrite	Urobilinogen	Protein
10 Leu/µL	0.05 mg/dL	0.2 mg/dL	9.0 mg/dL
Blood	Ketones	Bilirubin	Glucose
5 Ery/µL	2.5 mg/dL	0.4 mg/dL	60 mg/dL

Analytical Specificity				
Ammonium chloride	Ascorbic acid	Calcium chloride	Citric acid	
200 mg/dL	200 mg/dL	80 mg/dL	65 mg/dL	
Sodium nitrite	Oxalic acid	Sodium nitrate	Urea	
10 mg/dL	70 mg/dL	0.3 mg/dL	4000 mg/dL	
Lactose	Fructose	Galactose	Glycine	Creatinine
1.0 mg/dL	1.2 mg/dL	0.5 mg/dL	450 mg/dL	300 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	60
Total Agreement	98.5%

PRECISION:

It was determined using samples at 2 decision levels, one run in duplicate per day for 20 days. The results obtained were in agreement with those expected.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Only use strips whose packages are perfectly sealed.
- The test reading should not exceed 5 minutes for urine and 7 minutes for serum.
- Keep the bottle always closed and with the desiccant.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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