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PLASMAS COAGULATION CONTROL /  
PLASMAS CONTROLE DE COAGULACIÓN  
Ref. 36.003.00

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### ENGLISH

### INTENDED USE

Kit intended for internal quality control in Prothrombin Time (PT) and Activated Partial Thromboplastin Time (APTT) assays. Use in *in vitro* diagnostic.

### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.

- After reconstituted, the product in use is stable for up to 6 hours as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

### WORKING PRINCIPLE

Human lyophilized plasma used to verify the performance of the determination of Prothrombin Time (PT) and Activated Partial Thromboplastin Time (aPTT).

### PREPARATION AND HANDLING OF THE PRODUCT

Reconstitute with 1.0 mL of purified water. Homogenize gently by inversion, avoiding foaming. Wait 30 minutes at room temperature until complete dissolution of the product.

### PRODUCT DESCRIPTION

<b>CONTROL 1</b>	Plasma derivative lyophilized containing hemostasis analytes, plus buffer and preservatives.		
<b>CONTROL 2</b>	Plasma derivative lyophilized containing hemostasis analytes, plus buffer and preservatives.		

The controls present negative results for the HBs antigen and for the HIV and HCV antibodies. However, they should be treated as potentially infectious.

### NECESSARY EQUIPMENT FOR TESTING

- Glass and/or automatic pipettes.
- Coagulometer.
- Prothrombin Kit (PT - CAT BT 36.001.00).
- Activated Partial Thromboplastin Kit (aPTT - CAT BT 36.002.00).

### TEST PROCEDURE

The controls should be used in the same way as patient samples, according to the procedure established for the Prothrombin (TP) and Activated Partial Thromboplastin Time (aPTT) assays, and in accordance with the equipment protocol.

### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix controls from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the product if it displays any signs in disagreement with the ones specified in its MSDS.
- Avoid leaving controls outside the specified storage conditions.

### WARNINGS AND PRECAUTIONS

- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.bioteecnica.ind.br](http://www.bioteecnica.ind.br) or calling +55 3214 4646.
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.

### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 3214 4646, your local distributor or sending an e-mail to [sac@biotecnicaltda.com.br](mailto:sac@biotecnicaltda.com.br)

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	Plasmas Controle de Coagulação	CONTROLE 1	1 x 1 mL
	Plasmas Coagulation Control	CONTROLE 2	1 x 1 mL

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- FERNÁNDEZ-ROMERO, J.M.; CASTRO, M.D. Luque de. Continuous flow system for the evaluation of the extrinsic coagulation pathway. *Talanta*, v. 43, p.1531-1537, 1996.
- CLSI. One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline-Second Edition; Approved Guideline. **CLSI document H47-A2**. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- WHO Expert Committee on Biological Standardization. Guidelines for thromboplastins and plasma used to control oral anticoagulant therapy. World Health Organization, 1999: 7-8 (WHO Technical Report Series, N° 889).
- WANG, Huimin et al. Study on hemostatic mechanism of fully soluble hemostatic fiber. *Blood Coagulation And Fibrinolysis*, v. 18, n. 6, p.555-558, 2007.
- VAN HERWEGEN, Femke et al. The bleeding child. Part II: Disorders of secondary hemostasis and fibrinolysis. *Eur J Pediatr*. V. 171, p. 207-214, 2012.
- Brasil. Ministério da Saúde. Secretaria de Atenção à Saúde. Departamento de Atenção Especializada e Temática. **Manual de Diagnóstico Laboratorial das Coagulopatias Hereditárias e Plaquetopatias**. Brasília: Ministério da Saúde, 2016. 140 p.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests** - vol. 2, 5 ed. Washington DC: AACC Press, 2000.

**VALORES DOS PLASMAS CONTROLE DE COAGULAÇÃO / PLASMAS COAGULATION CONTROL VALUES / VALORES DE LOS PLASMAS CONTROLE DE COAGULACIÓN I Y II**

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These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

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