Colinesterase

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ENGLISH

INTENDED USE

Kit intended to determine the activity of Cholinesterase in serum and plasma samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Butyrylthiocholine Kinetic

Cholinesterase hydrolyzes butyrylthiocholine to form butyrate and thiocholine. Thiocholine reduces hexacyanoferrate III (yellow) to hexacyanoferrate II (colorless). The decrease in absorbance is spectrophotometrically determined at 405 nm, being directly proportional to the enzymatic activity.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (heparin and EDTA)

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	7 days
	-20 °C	1 year

PRODUCT DESCRIPTION

R 1 Sodium pyrophosphate ≥ 20 mmol / L; Potassium ferricyanide ≥ 0.5 mmol / L and stabilizer.

R 2 S-butrylcholine ≥ 20 mmol / L and stabilizers.

Traceable to specific absorptivity of chromogen hexacyanoferrate III using calibrated pipettes and handheld spectrophotometer that provide absolute values.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the control below:

Pathological Control – Quantial **REF** 13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Sample
Purified water	20 µL	-
Sample	-	20 µL
R1	-	1.0 mL
Homogenize and incubate at 37 °C for 3 minutes		
R2	200 µL	200 µL

2. Homogenize and insert into the cuvette holder thermostated at 37 °C, zeroed with purified water at 405nm. Set the stopwatch.
3. After 90 seconds, measure the initial absorbance A0 and read it again after 30, 60 and 90 seconds (A1, A2 and A3, respectively).

B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):

$$\Delta A/\text{min} = \frac{(A_0 - A_1) + (A_1 - A_2) + (A_2 - A_3)}{3} \times 2$$

The Cholinesterase activity is calculated by multiplying with the following equation:

Cholinesterase (U/L) = [(ΔA/min of Sample - ΔA/min of White)] x 65800

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

C) INTERPRETATION

Cholinesterase is an enzyme found in the liver, pancreas, heart, serum and the white matter of the brain, but its physiological role is unknown. Measurement of cholinesterase in serum is useful as an indicator of possible insecticide intoxication, for the detection of patients with atypical form of the enzyme, or as a liver function test. Decreased levels of cholinesterase are found in cases of organophosphorus poisoning, cirrhosis, myocardial infarction, acute infections and atypical phenotypes of the enzyme.

PERFORMANCE CHARACTERISTICS

Operating range
348.08 to 15000 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
5.71 U/L	348.08 U/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
1000 mg/dL	40 mg/dL	2000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 1.000x - 6.190$
Correlation Coefficient (R)	0.9944

By applying the regression equation, the total systematic error estimated for decision levels of 6000 U/L and 14000 U/L was -0.100% and -0.044%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (U/L)	Replicate s	Within-Run Precision		Total Precision	
		SD (U/L)	CV (%)	SD (U/L)	CV (%)
1520.40	80	7.300	0.50	27.790	1.80
4101.98	80	71.458	1.70	76.446	1.90
7163.69	80	37.637	0.50	45.814	0.60

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

37 °C	4900 - 11900 U/L
-------	------------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µKat/L):

Cholinesterase (U/L) x 0,0167 = Cholinesterase (µkat/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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