
DIRECTIONS FOR USE
Disposable Needle Electrodes

Applicable for:

- Disposable Concentric Needle Electrodes
- Disposable Hypodermic Needle Electrodes

Safe removal of protection sheath: twist and pull.

Technomed Europe
Amerikalaan 71
6199 AE Maastricht-Airport
The Netherlands

USA
Rx ONLY

STERILE EO

English

DISPOSABLE NEEDLE ELECTRODES

Intended use

Needle Electrodes are intended to be inserted in the subdermal, muscle or nerve tissue to sense bio-electric, EMG or EEG, signals distally and are intended to be proximal connected to electromyography / electroencephalogram recording equipment. The Disposable Hypodermic EMG Needle Electrode is used for muscle stimulation, motor unit action potential recording and drug delivery. The motor nerves are monitored by detecting EMG activity in the muscles they innervate. The drug used should be Botox Botulinum Toxin type A. Note: Technomed Europe does not supply any drugs with the Disposable Hypodermic EMG needle electrodes nor does Technomed Europe offer for sale any form of drugs.

Directions for use

When selecting a disposable needle electrode, it is important to use one that is of suitable length and diameter for the investigation being undertaken.

CAUTION!

- The disposable needle electrodes are sterile unless packaging is damaged or opened or use before date is passed.
- When introducing a disposable needle electrode into a patient it is recommended that it is not inserted up to the hub.
- If the disposable needle electrode bends before, during, or after insertion, do not straighten or re-insert. Bent disposable needle electrodes should be discarded and replaced with a new disposable needle electrode.
- Always discard used disposable needle electrodes in a properly marked medical biohazard sharps container.

WARNING!

These disposable needle electrodes are not intended for use as stimulating electrodes*) or as patient ground electrodes. If used as recording electrodes, the possible loss of electrocautery return path may result in high current density passing through these needles with possible tissue heating and damage resulting.

For single patient use only. Do not reuse or re-sterilize. Cleaning and re-sterilization of the device can affect the safety, performance and effectiveness of the device and expose patients and users to unnecessary risk.

*) The Disposable Hypodermic Needle Electrode can be used for muscle stimulation, recording EMG and drug delivery.

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Specifications:

Print : two sided black
Paper : 70 gr/m², white
Plano format: 17.0 x 45.5 cm
Finished : folded 6.5 x 17.0 cm (±0.2 cm)
Packed : manageable
Print: one front page (page 1) and 27 languages (page 2-28)

The diagram shows a table with two main sections. The left section is a 2x7 grid of cells. The top row of this section contains the numbers 1, 3, 5, 7, 9, 11, and 13. The bottom row contains the numbers 2, 4, 6, 8, 10, 12, and 14. The word "Front" is printed below the number 1. To the left of this grid is a vertical double-headed arrow labeled "17cm". Below the grid is a horizontal double-headed arrow labeled "6,5cm". To the right of the first grid is a second 2x7 grid of cells. The top row of this second grid contains the numbers 15, 17, 19, 21, 23, 25, and 27. The bottom row contains the numbers 16, 18, 20, 22, 24, 26, and 28.

1	3	5	7	9	11	13	15	17	19	21	23	25	27
2	4	6	8	10	12	14	16	18	20	22	24	26	28