PCR LÁTEX

CRP LATEX / PCR LÁTEX
Ref. 30.005.00 (KIT) / 30.005.00 (REAGENTE)

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ENGLISH

INTENDED USE

Kit for the determination of C-reactive Protein (CRP) in serum. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Ready-to-use reagents and controls.

WORKING PRINCIPLE

Method: Plate agglutination

Latex PCR is a plate agglutination method for qualitative and semi-quantitative detection, where the latex particles coated with anti-human PCR antibodies are bound by the PCR present in the sample

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

| | Temperature | Stability Period |
|-------|-------------|------------------|
| Serum | 2 to 8 °C | 7 days |
| | -20 °C | 3 months |

PRODUCT DESCRIPTION

Suspension of Latex particles coated with goat anti-human IgG, pH 8.2. Preservative. The sensitivity of RLAX is traceable to ERM-DA 472/IFCC (International Federation of Clinical Chemistry)

Human serum with a concentration of CRP> 20 mg/L. Preservative. (INCLUDED IN THE PCR LATEX KIT)

Human serum in buffer Phosphate, stabilizer and preservative. (INCLUDED IN THE PCR LATEX KIT)

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the controls below:

| | | |
|---|------------|-----------|
| CONTROL (+) (INCLUDED IN THE PCR LATEX KIT) | REF | 30.005.00 |
| CONTROL (-) (INCLUDED IN THE PCR LATEX KIT) | | 30.005.00 |

NECESSARY EQUIPMENT FOR TESTING

- Reaction plate and disposable rod;
- Pipette, clock or stopwatch..

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Automation: This procedure **cannot** be applied to automated analyzers.

Qualitative Method

- Before starting the tests, allow the RLAX, controls and samples to reach room temperature. Homogenize before use.
- On the reaction plate use a different circle for each sample or control.
- Place 20 µL of the sample / controls.
- In each circle add 20 µL of RLAX.
- Mix each determination using disposable rods (one for each sample).
- Homogenize the plate for 2 minutes in circular motions, tilting it back and forth. Observe the result under artificial light: Presence or absence of agglutination.

Semi-quantitative method

- Dilute the sample with saline (0.9% NaCl) as follows:

| Dilution | Saline (µL) | Sample(µL) | Dilution Factor |
|----------|-------------|---------------|-----------------|
| 1/2 | 100 | 100 (pure) | 2 |
| 1/4 | 100 | 100 (sol.1/2) | 4 |
| 1/8 | 100 | 100 (sol.1/4) | 8 |
| 1/16 | 100 | 100 (sol.1/8) | 16 |

- Carry out the test of each dilution according to the qualitative method.

B) CALCULATIONS

The approximate concentration in the semi-quantitative test is obtained by multiplying the dilution factor by the sensitivity of the test (6 mg/L).

Exemplo:

Agglutination up to 1/8 dilution.

Dilution Factor = 8

PCR = 48 IU / mL (8 x 6)

C) INTERPRETATION

C-Reactive Protein is an acute phase protein that can significantly increase in most bacterial and viral infectious processes, inflammation and malignant neoplasms. Its elevation occurs rapidly and can reach values of 300 mg / L in 12 to 24 hours.

PERFORMANCE CHARACTERISTICS

| Operating range |
|-----------------|
| Up to 300 mg/L |

| Sensitivity |
|-----------------|
| Detection Limit |
| 6 mg/L |

| Analytical Specificity | | | |
|------------------------|-----------|---------------|-------------------|
| Hemoglobin | Bilirubin | Triglycerides | Rheumatoid Factor |
| 1000 mg/dL | 40 mg/dL | 2000 mg/dL | 40 UI/mL |

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

| Accuracy | |
|--------------------|------|
| Number of Samples | 100 |
| Positive Agreement | 96% |
| Negative Agreement | 100% |

Prozone High Dose Effect:

There is no prozone effect up to 300 mg/L of PCR.

Precision:

Intra-assay: the study was carried out in two daily runs, in duplicate, for 10 days, using 04 samples, two positive and two negatives, with 100% agreement.

Inter-assay: the study was carried out in a daily run of four replicates, for 5 days, using 04 samples, two positive and two negatives, with 100% agreement

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

| | |
|--------|--------------|
| Adults | Up to 6 mg/L |
|--------|--------------|

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicaltda.com.br.

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- BURTIS, C.A.; ASHWOOD, E.R.; BRUNS, D.E. **Tietz:** Fundamentos de Química Clínica. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- HOKAMA, Yoshitsugi; NAKAMURA, Robert M. C-Reactive Protein: Current Status and Future Perspectives. *Journal Of Clinical Laboratory Analysis*, v. 1, p.15-27, 1987.
- OTSUJI, Shogo; SHIBATA, Hideaki; UMEMA, Mamoru. Turbidimetric Immunoassay of Serum C-reactive Protein. *Clin. Chem.* v. 28, p.2121-2124, 1982.
- PEPYS, M. B.; DASH, A. C.; ASHLEY, M. J. Isolation of C-reactive protein by affinity chromatography. *Clin. Exp. Immunol.* v. 30, p.32-37, 1977.
- WADSWORTH, Charles; WADSWORTH, Elsa. Efficacy of latex agglutination and quantification methods for determination of C-reactive protein (CRP) in pediatric sera. *Clinica Chimica Acta*, v. 138, p.309-318, 1984.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- YAMAMOTO, S. et al. Avidity of antibody and agglutinability of antibody-sensitized latex in latex agglutination test. *Veterinary Immunology And Immunopathology*, v. 36, p.257-264, 1993.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: **AAAC Press**, 2000.

| TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES | | | |
|--|--|--------------|--|
| | Consultar as instruções para utilização Consult instructions for use Consulté las instrucciones de uso | | Descartar corretamente Dispose properly Desechar adecuadamente |
| REF | Número de catálogo Catalog number Número de catálogo | R <N> | Reagente Reagent Reactivo |
| LOT | Código do lote Batch code Código de lote | | Límite de temperatura Temperature limitation Límite de temperatura |
| | Risco biológico Biological risk Riesgo biológico | X | Nocivo / Irritante Harmful / Irritant Nocivo / Irritante |
| CONTROL | Controle Control Control | | Validade Use by date Fecha de Caducidad |
| IVD | Produto para a saúde para diagnóstico <i>in vitro</i> In Vitro Diagnostic medical device Producto sanitario para diagnóstico <i>in vitro</i> | | |